
Operating Manual

Body Fat Scale

FAT 200

BODY FAT MONITOR, OPERATION

FAT200

1. WEIGHING MODE



Turn the scale on by pressing the  button. After the initial **HELO** message, the display flashes **1234**. During this period press the  key a second time to enter weighing mode. The display will read **0.0**. The scale can be set to read in kg (kilograms) or lb (pounds) by pressing the  button for one second and releasing it.

The scale will turn off automatically after a few minutes.

2. PROGRAMMING PERSONAL DATA.

Turn the scale on with the  button. When the display flashes **1234**, press the  button and choose one of the eight memory storage locations by use of the  and  buttons followed by the  button.


The display now shows a male or female symbol; choose one by use of the  and  buttons followed by the  button.

The display now shows height in cm (centimeters). Enter your personal height by use of the  and  buttons followed by the  button. If you do not know your height in cm, multiply your height in inches by 2.54.

The display now shows age in years. Enter the nearest age by use of the  and  buttons followed by the  button.

Personal data entry is complete when the display flashes **1234** again. At this time you can repeat the process by pressing  again or enter the body fat weighing mode by pressing  or .

3. MEASURING BODY FAT

Turn the scale on with the  button. When the display flashes **1234**, press the  or  button until the display shows your personal memory number and wait until the display shows **0.0**. Change the units to kg or lb by pressing the  button for one second and releasing it. Remove shoes and socks, then step on the scale and wait. The display first shows your weight, followed by flashing **0000** during the body resistance measurement period. This is followed by an alternating display of your weight and your body fat percentage.

4. BATTERY INSERTION/REPLACEMENT

Snap open the battery cover on the underside of the scale and remove the old batteries. Replace with two AA size batteries, taking care to insert them with the + and – as shown, - end towards the spring.

WARNING : People with pacemakers or other electronic body implants should not use this product since the body resistance measurement method may interfere with its proper operation.

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